

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

SURFACE DISINFECTION AND CLEANING OF INVASIVE MEDICAL DEVICES WITH MECHANICAL ACTION

1 / 9

Broad spectrum* of activity with short contact times**

High material and product compatibility

due to oxidative base

Suitable for use in the food sector



HACCP

* bactericidal, levurocidal, tuberculocidal, mycobactericidal, sporicidal, fungicidal and virucidal

** Surface disinfection, tuberculocidal surface disinfection, disinfection against Clostridium difficile within 5 minutes, sporicidal surface disinfection within 15 minutes.

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

PRODUCT DESCRIPTION

2 / 9

ULTRASOL OXY WIPES are wipes for surface disinfection and cleaning of invasive medical devices with mechanical action. The wipes have a broad spectrum of activity* against bacteria and viruses, including spores.

Thanks to its oxidative base, **ULTRASOL OXY WIPES** have a high level of material compatibility and can be used on almost all materials. It leaves no toxic or environmentally harmful residues on the surface.

Our **ULTRASOL OXY WIPES XL** are suitable for disinfecting large surfaces due to the larger wipes. The pre-saturated wipes allow convenient and time-saving application, even when used in difficult conditions.

ULTRASOL OXY WIPES are effective against bacteria, yeasts (including *Candida auris*) and enveloped viruses.

APPLICATIONS AND NOTES

According to Biocidal Products Regulation (BPR)

For rapid disinfection and cleaning of alcohol sensitive medical equipment and surfaces of every type. Suitable for use in the food sector.

Acc. to EU Medical Devices Regulation

Surface disinfection and cleaning of invasive medical devices with mechanical action.

Further areas of application

In addition to the medical sector also suitable for the food sector and large canteen kitchens as well as for industry and public facilities.

Application

For professional use only. Visible soiling must be removed before disinfection. Wipe the surfaces to be disinfected thoroughly. Ensure that the surface is completely wet throughout the entire exposure time.

When disinfecting, the exposure time before reuse must be taken into account. Use personal protective equipment (protective gloves). Avoid direct skin contact. According to the EU Medical Device Regulation, users/patients are obligated to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user/patient is established.

Application notes

Shelf life after opening: ULTRASOL OXY WIPES: 42 days
ULTRASOL OXY WIPES XL: 42 days

When using, please follow the instructions provided by the medical device manufacturer. When disinfecting incubators for premature infants, the KRINKO guidelines must be observed.

Composition

100 g solution contain: 7 g Hydrogen Peroxide, 0.1 g Peracetic Acid, 0.1 g Glycolic Acid.

Material compatibility

Wide range of applications on surfaces and medical devices.
(see page 4 - 7)

Product status

Dual Registration (medical device/biocide)

Precautionary and hazard statements

Causes serious eye irritation. Wear protective gloves. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Dispose of contents/container to approved disposal company or local collection.

Use disinfectants safely.

Always read label and product information before use.

* bactericidal, levurocidal, tuberculocidal, mycobactericidal, sporicidal, fungicidal and virucidal

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

SPECTRUM OF EFFICACY AND CONTACT TIMES

3 / 9



SPECTRUM OF ACTIVITY AND CONTACT TIMES *			1 min	5 min	15 min
recommendation for surface disinfection					
bactericidal, levurocidal	EN 13624 EN 13727 EN 16615	clean and dirty conditions		•	
tuberculocidal	EN 14348	clean and dirty conditions		•	
mycobactericidal	EN 14348	clean and dirty conditions		•	
sporicidal against C. diff. ULTRASOL OXY WIPES	EN 17126 EN 17846	dirty conditions		•	
sporicidal against C. diff. ULTRASOL OXY WIPES XL	EN 17126 EN 17846	dirty conditions		•	
sporicidal ULTRASOL OXY WIPES	EN 17126 EN 17846	clean and dirty conditions			•
sporicidal ULTRASOL OXY WIPES XL	EN 17126 EN 17846	dirty conditions			•
fungicidal	EN 13624	clean and dirty conditions		•	
virucidal	EN 14476	clean and dirty conditions		•	
limited virucidal	EN 14476	clean and dirty conditions		•	
additional test results					
bactericidal	EN 1276	dirty conditions		•	
yeastcidal (Candida albicans)	EN 1650	dirty conditions		•	
sporicidal	EN 13704	dirty conditions			•

* The spectrum of activity and contact times apply both to use as a biocide and as a medical device.

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

MATERIAL COMPATIBILITY

4 / 9



MATERIAL METALS	not recommended	recommended	APPLICATION PRODUCT EXAMPLE
stainless steel V2A		<ul style="list-style-type: none"> Rolling walking aids Toilet chairs 	
aluminum anodized	<ul style="list-style-type: none"> 		<ul style="list-style-type: none"> Walking frames

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

MATERIAL COMPATIBILITY

5 / 9



MATERIAL PLASTICS: ELASTOMERS	not recommended	recommended	APPLICATION PRODUCT EXAMPLE
silicones		<ul style="list-style-type: none"> • 	<ul style="list-style-type: none"> • Face masks
PUR (polyurethane)		<ul style="list-style-type: none"> • 	<ul style="list-style-type: none"> • Open cuff face mask
CR (neoprene)		<ul style="list-style-type: none"> • 	<ul style="list-style-type: none"> • Resuscitator bag
EPDM (ethylene propylene diene (monomer) rubber)		<ul style="list-style-type: none"> • 	<ul style="list-style-type: none"> • Nursing trolleys
TPS (styrene TPE)		<ul style="list-style-type: none"> • 	
NBR (nitrile butadiene rubber)		<ul style="list-style-type: none"> • 	

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

MATERIAL COMPATIBILITY

6 / 9



MATERIAL PLASTICS: THERMOPLASTICS	not recommended	recommended	APPLICATION PRODUCT EXAMPLE
PC (polycarbonate e.g. Makrolon)		•	Ultrasound devices
		•	MRI devices
		•	EEG devices
		•	ECG devices
		•	CT devices
PC/ABS (polycarbonate/acrylonitril-butadiene-styrene)		•	X-ray devices
		•	Ultrasound probes e.g. transvaginal and abdominal probes
		•	Incubators
ABS (acrylonitril-butadiene-styrene)		•	Patient monitoring monitors
PEI (polyetherimide)		•	Sterilization and transport containers
PMMA (polymethylmethacrylate)		•	Acrylic and plexiglass incubators
PA (polyamide)		•	X-ray apron
PE-HD (polyethylene-high density)		•	Storage and transport containers
PP (polypropylene)		•	Hose assemblies
PVC (polyvinylchloride)		•	Oxygen bag
		•	Bag for training manikin
		•	Emergency bag

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY[®] WIPES



PROTECT

MATERIAL COMPATIBILITY

7 / 9



MATERIAL PLASTICS: TEXTILES	not recommended	recommended	APPLICATION PRODUCT EXAMPLE
Lyocell	•		
Viscose	•		
Cotton	•		
textile fabric	•		

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY[®] WIPES



PROTECT

PURCHASING INFORMATION

8 / 9

Product	Single unit	Unit	Content	Wipe size	REF
ULTRASOL OXY WIPES	Package	6	72 Wipes	20 x 20 cm	00-270-To72
ULTRASOL OXY WIPES XL	Dispenser System	4	120 Wipes	17,5 x 36 cm	00-270-OSEB120

National information may differ. For further information, please contact our subsidiary or your local dealer.



CERTIFICATIONS



For an overview of our current certifications, scan the QR code to visit:

<https://www.schumacher-online.com/en/cert>

Dr. Schumacher is certified according to DIN EN 13485, DIN EN ISO 9001, DIN EN ISO 14001, BS OHSAS 18001, has a validated environment management system according to EMAS and is a member of IHO, VCI, BAH, DGSV and of the DGKH.

SURFACE / READY TO USE DISINFECTANTS

ULTRASOL OXY® WIPES



PROTECT

PRODUCT FAMILY OVERVIEW

9 / 9



ULTRASOL OXY® WIPES



ULTRASOL OXY® WIPES XL



ULTRASOL OXY®



RECOMMENDED NON-WOVEN WIPE DISPENSER SYSTEMS



ONE SYSTEM+ PLUS



ONE SYSTEM BASIC



DESCO WIPES



ECO WIPES