Dr.
Schumacher

Wir schützen
Deine Gesundheit.

Dosiertabelle / Dosage table

Fertige Gebrauchslösung Ready-to-use solution	0,25 %	0,5 %	1,0 %	1,5 %	2,0 %	3,0 %	4,0 %
1 L	2,5 ml	5 ml	10 ml	15 ml	20 ml	30 ml	40 ml
2 L	5 ml	10 ml	20 ml	30 ml	40 ml	60 ml	80 ml
3 L	7,5 ml	15 ml	30 ml	45 ml	60 ml	90 ml	120 ml
4 L	10 ml	20 ml	40 ml	60 ml	80 ml	120 ml	160 ml
5 L	12,5 ml	25 ml	50 ml	75 ml	100 ml	150 ml	200 ml
8 L	20 ml	40 ml	80 ml	120 ml	160 ml	240 ml	320 ml
10 L	25 ml	50 ml	100 ml	150 ml	200 ml	300 ml	400 ml
16 L	40 ml	80 ml	160 ml	240 ml	320 ml	480 ml	640 ml
20 L	50 ml	100 ml	200 ml	300 ml	400 ml	600 ml	800 ml
30 L	75 ml	150 ml	300 ml	450 ml	600 ml	900 ml	1200 ml

Herstellung gebrauchsfertiger Desinfektionslösungen:
Die Anwendungslösung ist aus Konzentrat und Leitungswasser herzustellen. Um 5 Liter Anwendungslösung mit einer 2%igen Konzentration herzustellen, müssen 100 ml des Konzentrats zu 4,9 Liter Leitungswasser hinzugefügt werden. Die Richtlinien und Einwirkzeiten des Herstellers sind zu beachten.

Preparation of a ready-to-use solution of disinfectants:
The ready-to-use solution consists of concentrate and tap water. In order to prepare 5 liters ready-to-use solution with a concentration of 2 % you have to add 100 ml of concentrate to 4,9 liters of tap water. Manufacturers guidelines and exposure times should be followed.

Dr. Schumacher GmbH • Am Roggenfeld 3 • 34323 Malsfeld • T +49 5664 9496-0 • info@schumacher-online.com

www.schumacher-online.com

© Dr. Schumacher GmbH · PSP-000611 · Stand 03/2026